



How Do I Use the Entyvio® (vedolizumab) Pen?

Before use:




After use:



1. Get everything ready

- Find a relatively calm spot – you'll want to give yourself some time and space to get things right
- Take one injection pen box from the fridge so it can start to warm up to room temperature (this takes about 30 minutes)

- Grab the other supplies you need:
 - An alcohol pad
 - Cotton wool ball or gauze
 - Sharps disposal container



2. Check the injection pen

- Wash your hands
- Look over the injection pen for any broken or missing seals
- Check the expiry date
- Check the medicine looks clear and colourless, or light yellow

- Don't worry if there are air bubbles in the pen – this is normal, don't shake it



3. Pick an injection site





- Choose** where you will inject – front of thighs, tummy (not around the belly button) or back of the upper arm (if someone else is doing the injection)
- Make sure you **inject into a different place** from last time you injected, it may be useful to keep a diary

- Wipe the area** you plan to inject with an **alcohol pad** (don't touch this area again before you inject)




4. Injecting Entyvio®

- Pull the purple cap off the injection pen**
- Hold the pen so you can see the viewing window
- Place the yellow end of the pen at 90-degrees** to the area you cleaned earlier (it should be straight up and down)
- Push down on the pen as far as it will go** to begin the injection
- Hold and count to 10** while pushing down as firmly as you can. This will allow all the medicine to be injected
- You might hear two clicks, one at the start and one near the end
- Check that the viewing window is filled with purple** before you stop pushing
 - It's fine if you see a small bit of grey in the window
 - Lift the pen from the injection site
- The yellow needle shield will drop down and lock over the needle
- If the viewing window isn't completely full, you might not have received your full dose – contact your healthcare team if this happens
- If you see a small amount of blood at the injection site, press on your skin with the cotton wool ball or gauze



5. Clean up after yourself

- Put the **injection pen into a sharps container** straight away – don't put the syringe in your normal rubbish. You should have received a sharps container with your first delivery of Entyvio® subcutaneous (SC)
- If you don't have one of those, please contact your homecare provider* immediately and use a container that:
 - Is made of heavy-duty plastic
 - Can be closed with a tight-fitting, puncture-resistant lid without anything being able to come out
 - Is upright and stable during use
- Is leak-resistant
- Is properly labelled to warn of hazardous waste inside the container
- You can put the rest of the things (the cap, the alcohol wipe, the cotton wool ball or gauze) in your normal rubbish or in the sharps bin
- When your sharps bin is almost full, please contact your homecare provider* who will remove your sharps disposal container and provide you with a new one



Please read the patient information leaflet, which you can find in your Entyvio® packaging, for further information about your treatment.